

# LABORATORY / INDUSTRIAL PLASTIC BODY CONDUCTIVITY CELL TCP 01 & TCP 1



### SPECIFICATIONS

	TCP 01	TCP 1
1. Electrode Material	: Graphite	Graphite
2. No. of Electrode	: 2 Electrodes	2 Electrodes
3. Body Material	: Plastic	Plastic
4. Cell Constant	: 0.1 cm <sup>-1</sup>	1 cm <sup>-1</sup>
5. Measuring Range	: 0-2 mS	0-20 mS
6. Sensor Diameter	: 12 mm	12 mm
7. Temp. Range	: -5 to 80°C	-5 to 80°C
8. Temp. Compensation	: Pt100/ Pt1000/ 30K NTC (Optional)	Pt100/ Pt1000/ 30K NTC (Optional)
9. Connector	: BNC & Lugs	BNC & Lugs

#### **OTHER INDUSTRY:**

Pharma, Petrochemical, Swimming Pools, Aquarium, Electroplating, Food, Chemical, Fertilizer, Dairy, Oil, Drinking Water, Ground Water, Waste Water

#### **OTHER APPLICATIONS:**

Extreme pH (Acidic), Extreme pH (Alkali), Narrow Vessel/ Semi Micro, Low Ionic/ Ultrapure, Oil/ Water Emulsion, HF, Low Temp. (-30°C), High Temp. (135°C), Steam Sterilisable

**TOSHNIWAL INSTRUMENTS MANUFACTURING PVT. LTD.,**

P.O. Gagwana - 305023, Distt. Ajmer (Raj.) INDIA

Phone: +91-145-2971131/2/3; Mobile: +91-9829112345

Email: [info@toshcon.com](mailto:info@toshcon.com), [support@toshcon.com](mailto:support@toshcon.com); Website: [www.toshcon.com](http://www.toshcon.com)

REGD OFFICE: 401, Manish Chamber, Sonawala Cross Lane, Goregaon (E), MUMBAI - 400 063